

CAB19187

Product Information

Product SKU: CAB19187	Gene ID: 107276289	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Oryza sativa

Additional Information

Observed MW: 50kDa	Conjugate: -
Calculated MW: 39kDa	Isotype: IgG

Immunogen Information

Background: Transcriptional activator that acts as floral inducer to promote short-day (SD) flowering pathway. Activates HD3A and other FT-like genes independently from HD1. May also activate MADS-box transcription factors involved in flowering regulation . Functions as a response regulator involved in His-to-Asp phosphorelay signal transduction system. Phosphorylation of the Asp residue in the receiver domain activates the ability of the protein to promote the transcription of target genes. May directly activate some type-A response regulators in response to cytokinins

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: ef1; Ehd1; RR30; OsJ_31808; EHD1

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa EHD1 .

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.