

CAB19250

Product Information

Product SKU:	CAB19250	Gene ID:	1796	Size:	20uL, 100uL
Clone No:	ARC2410	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	62kDa	Conjugate:	Unmodified
Calculated MW:	52kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is part of a signal transduction pathway downstream of receptor tyrosine kinases. The encoded protein is a scaffold protein that helps form a platform for the assembly of multiprotein signaling complexes. Several transcript variants encoding different isoforms have been found for this gene.

Recommended Dilution: WB,1:500 - 1:1000 IF/ICC,1:50 - 1:200

Synonyms: pp62; P62DOK; DOK1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 382-481 of human DOK1 (Q99704).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.