

CAB19260

Product Information

Product SKU:	CAB19260	Gene ID:	26052	Size:	20uL, 100uL
Clone No:	ARC2428	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	98kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of a family of guanosine triphosphate (GTP)-binding proteins that associate with microtubules and are involved in vesicular transport. The encoded protein functions in the development of megakaryocytes. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: Dyna III; Dynamin 3

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 700-800 of human Dynamin 3 (Q9UQ16).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.