

CAB19261

Product Information

Product SKU:	CAB19261	Gene ID:	5309	Size:	20uL, 100uL
Clone No:	ARC2429	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	32kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the RIEG/PITX homeobox family, which is in the bicoid class of homeodomain proteins. Members of this family act as transcription factors. This protein is involved in lens formation during eye development. Mutations of this gene have been associated with anterior segment mesenchymal dysgenesis and congenital cataracts.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: ASMD; ASOD; PTX3; ASGD1; CTPP4; CTRCT11; PITX3

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PITX3 (O75364).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.