

## CAB19263

---

### Product Information

<b>Product SKU:</b>	CAB19263	<b>Gene ID:</b>	873	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2432	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse

---

### Additional Information

<b>Observed MW:</b>	33kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	30kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** The protein encoded by this gene belongs to the short-chain dehydrogenases/reductases (SDR) family, which function as NADPH-dependent oxidoreductases having wide specificity for carbonyl compounds, such as quinones, prostaglandins, and various xenobiotics. Alternatively spliced transcript variants have been found for this gene.

**Recommended Dilution:** WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

**Synonyms:** CBR; hCBR1; PG-9-KR; SDR21C1; CBR1

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 178-277 of human CBR1 (P16152).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.