

CAB19265

Product Information

Product SKU:	CAB19265	Gene ID:	9355	Size:	20uL, 100uL
Clone No:	ARC2435	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	44kDa	Conjugate:	Unmodified
Calculated MW:	44kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a protein belonging to a large protein family, members of which carry the LIM domain, a unique cysteine-rich zinc-binding domain. The encoded protein may function as a transcriptional regulator. The protein can recapitulate or rescue phenotypes in Drosophila caused by a related protein, suggesting conservation of function during evolution.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: LH2; hLhx2; LHX2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LHX2 (P50458).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.