

CAB19273

Product Information

Product SKU:	CAB19273	Gene ID:	55824	Size:	20uL, 100uL
Clone No:	ARC2448	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	76kDa	Conjugate:	Unconjugated
Calculated MW:	47kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a type III transmembrane adaptor protein that binds to the tyrosine kinase csk protein. It is thought to be involved in the regulation of T cell activation.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: CBP; PAG; PAG1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 333-432 of human PAG1 (Q9NWX8).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.