

CAB19277

Product Information

Product SKU:	CAB19277	Gene ID:	2516	Size:	20uL, 100uL
Clone No:	ARC2457	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	48kDa	Conjugate:	Unmodified
Calculated MW:	52kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene is a transcriptional activator involved in sex determination. The encoded protein binds DNA as a monomer. Defects in this gene are a cause of XY sex reversal with or without adrenal failure as well as adrenocortical insufficiency without ovarian defect.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: ELP; SF1; FTZ1; POF7; SF-1; AD4BP; FTZF1; SPGF8; SRXX4; SRXY3; hSF-1; Steroidogenic factor-1 (SF-1/NR5A1)

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 351-461 of human Steroidogenic factor-1 (SF-1/NR5A1) (Q13285).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.