

CAB19284

Product Information

Product SKU:	CAB19284	Gene ID:	3859	Size:	20uL, 100uL
Clone No:	ARC2465	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	54kDa	Conjugate:	Unconjugated
Calculated MW:	54kDa	Isotype:	IgG

Immunogen Information

Background: KRT12 encodes the type I intermediate filament chain keratin 12, expressed in corneal epithelia. Mutations in this gene lead to Meesmann corneal dystrophy.

Recommended Dilution: WB,1:100 - 1:500 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: K12; MECD1; Cytokeratin 12 (KRT12)

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 250-494 of human Cytokeratin 12 (KRT12) (NP_000214.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.