

CAB19287

Product Information

Product SKU:	CAB19287	Gene ID:	25898	Size:	20uL, 100uL
Clone No:	ARC2469	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	30kDa	Conjugate:	Unmodified
Calculated MW:	30kDa	Isotype:	IgG

Immunogen Information

Background: The protein encoded by this gene has ubiquitin ligase activity. It mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including tumor protein 53, histone deacetylase 1, and cyclin-dependent kinase inhibitor 1B, thus regulating their levels and cell cycle progression. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: ZCHY; ARNIP; CHIMP; PIRH2; RNF199; ZNF363; PRO1996; RCHY1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 160-261 of human RCHY1 (Q96PM5).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.