

CAB19292

Product Information

| | | | | | |
|---------------------|----------|----------------------|--------|--------------------|-----------------|
| Product SKU: | CAB19292 | Gene ID: | 79370 | Size: | 20uL, 100uL |
| Clone No: | ARC2474 | Host Species: | Rabbit | Reactivity: | Human,Mouse,Rat |

Additional Information

| | | | |
|-----------------------|-------|-------------------|--------------|
| Observed MW: | 37kDa | Conjugate: | Unconjugated |
| Calculated MW: | 37kDa | Isotype: | IgG |

Immunogen Information

Background: The protein encoded by this gene belongs to the BCL2 protein family. BCL2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. Overexpression of this gene has been shown to induce apoptosis in cells. Three alternatively spliced transcript variants encoding two distinct isoforms have been reported for this gene.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: BCLG; BCL2L14/BCLG

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human BCL2L14/BCLG (Q9BZR8).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.