

CAB19331

Product Information

Product SKU: CAB19331	Gene ID: 5458	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human, Mouse

Additional Information

Observed MW: Refer to figures	Conjugate: Unconjugated
Calculated MW: 43kDa	Isotype: IgG

Immunogen Information

Background: The protein encoded by this gene is a member of the POU-domain transcription factor family and may be involved in maintaining visual system neurons in the retina. The level of the encoded protein is also elevated in a majority of breast cancers, resulting in accelerated tumor growth.

Recommended Dilution: WB, 1:500 - 1:1000 IF/ICC, 1:50 - 1:200

Synonyms: BRN3B; BRN3.2; Brn-3b; POU4F2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human POU4F2 (NP_004566.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.