

## CAB19467

---

### Product Information

<b>Product SKU:</b>	CAB19467	<b>Gene ID:</b>	22828	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	170kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	141kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Enables RNA binding activity and RNA polymerase II C-terminal domain phosphoserine binding activity. Involved in negative regulation of termination of RNA polymerase II transcription, poly(A)-coupled and positive regulation of DNA-templated transcription, elongation. Located in nucleoplasm.

**Recommended Dilution:** WB,1:500 - 1:1000

**Synonyms:** RBM16; SCAF8

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 572-634 of human SCAF8 (NP\_055707.3).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.