

CAB19660

Product Information

Product SKU:	CAB19660	Gene ID:	1506	Size:	20uL, 100uL
Clone No:	ARC2218	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	28kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a serine-type endopeptidase with chymotrypsin- and elastase-2-like activities. The gene encoding this zymogen is expressed specifically in the pancreas and likely functions as a digestive enzyme.

Recommended Dilution: WB,1:5000 - 1:8000 IHC-P,1:50 - 1:200 IF/ICC,1:50 - 1:200

Synonyms: CTRL1; CTRL

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 165-264 of human CTRL (P40313).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.