

## CAB19711

---

### Product Information

<b>Product SKU:</b>	CAB19711	<b>Gene ID:</b>	3500	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2241	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	50kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	52kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to enable antigen binding activity and immunoglobulin receptor binding activity. Predicted to be involved in several processes, including activation of immune response; defense response to other organism; and phagocytosis. Predicted to act upstream of or within several processes, including immunoglobulin mediated immune response; positive regulation of hypersensitivity; and positive regulation of phagocytosis. Located in extracellular space.

**Recommended Dilution:** WB,1:1000 - 1:5000 IHC-P,1:100 - 1:500

**Synonyms:** -

**Purification Method:** Affinity purification

**Immunogen:** Human IgG

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.