

## CAB19713

---

### Product Information

<b>Product SKU:</b>	CAB19713	<b>Gene ID:</b>	3502	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2243	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human

---

### Additional Information

<b>Observed MW:</b>	60kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	41kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** Predicted to enable antigen binding activity and immunoglobulin receptor binding activity. Involved in retina homeostasis. Located in blood microparticle and extracellular exosome.

**Recommended Dilution:** WB, 1:1000 - 1:5000

**Synonyms:** IgG3; Human IgG3

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 50-150 of human IgG3 (P01860).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.