

CAB19720

Product Information

Product SKU:	CAB19720	Gene ID:	3514	Size:	20uL, 100uL
Clone No:	ARC2247	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	25kDa	Conjugate:	Unconjugated
Calculated MW:	12kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable antigen binding activity and immunoglobulin receptor binding activity. Involved in retina homeostasis. Located in blood microparticle and extracellular exosome.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: Km; HCAK1; IGKCD; Human Ig Kappa Chain

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-107 of human Ig Kappa Chain (P01834).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.