

CAB19734

Product Information

Product SKU:	CAB19734	Gene ID:	4824	Size:	20uL, 100uL
Clone No:	ARC2261	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	26kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a homeobox-containing transcription factor. This transcription factor functions as a negative regulator of epithelial cell growth in prostate tissue. Aberrant expression of this gene is associated with prostate tumor progression. Alternate splicing results in multiple transcript variants of this gene.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: NKX3; BAPX2; NKX3A; NKX3.1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-234 of human NKX3.1 (Q99801).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.