

CAB19736

Product Information

Product SKU:	CAB19736	Gene ID:	5004/5005	Size:	20uL, 100uL
Clone No:	ARC2263	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	45kDa	Conjugate:	Unconjugated
Calculated MW:	47kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a key acute phase plasma protein. Because of its increase due to acute inflammation, this protein is classified as an acute-phase reactant. The specific function of this protein has not yet been determined; however, it may be involved in aspects of immunosuppression.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: AGP-B; AGP-B'; AGP2; ORM1/2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human ORM1/2 (NP_000598.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.