

CAB19737

Product Information

Product SKU:	CAB19737	Gene ID:	5013	Size:	20uL, 100uL
Clone No:	ARC2264	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	37kDa	Conjugate:	Unconjugated
Calculated MW:	37kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a member of the bicoid sub-family of homeodomain-containing transcription factors. The encoded protein acts as a transcription factor and may play a role in brain and sensory organ development. A similar protein in mouse is required for proper brain and sensory organ development and can cause epilepsy. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: -

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human OTX1 (NP_055377.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.