

CAB19766

Product Information

Product SKU:	CAB19766	Gene ID:	10926	Size:	20uL, 100uL
Clone No:	ARC2293	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	77kDa	Conjugate:	Unconjugated
Calculated MW:	77kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable protein serine/threonine kinase activator activity. Predicted to be involved in positive regulation of nuclear cell cycle DNA replication and regulation of cell cycle phase transition. Located in nuclear body.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: ASK; CHIF; DBF4A; ZDBF1; DBF4

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human DBF4 (Q9UBU7).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.