

## CAB19768

---

### Product Information

<b>Product SKU:</b>	CAB19768	<b>Gene ID:</b>	22924	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	ARC2305	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

### Additional Information

<b>Observed MW:</b>	32kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	32kDa	<b>Isotype:</b>	IgG

---

### Immunogen Information

**Background:** The protein encoded by this gene is a member of the RP/EB family of genes. The protein localizes to the cytoplasmic microtubule network and binds APCL, a homolog of the adenomatous polyposis coli tumor suppressor gene.

**Recommended Dilution:** WB,1:500 - 1:2000 IF/ICC,1:50 - 1:200

**Synonyms:** EB3; RP3; EBF3; EBF3-S; EB3/MAPRE3

**Purification Method:** Affinity purification

**Immunogen:** A synthetic peptide corresponding to a sequence within amino acids 182-281 of human EB3/MAPRE3 (Q9UPY8).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.