

CAB19783

Product Information

Product SKU:	CAB19783	Gene ID:	23369	Size:	20uL, 100uL
Clone No:	ARC2308	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	110-125kDa	Conjugate:	Unconjugated
Calculated MW:	114kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a protein that belongs to a family of RNA-binding proteins. The encoded protein functions as a translational repressor during embryonic development and cell differentiation. This protein is also thought to be a positive regulator of cell proliferation in adipose-derived stem cells. Alternate splicing results in multiple transcript variants.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: PUMH2; PUML2; Pumilio 2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Pumilio 2 (Q8TB72).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.