

CAB19816

Product Information

Product SKU:	CAB19816	Gene ID:	140609	Size:	20uL, 100uL
Clone No:	ARC2342	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	32kDa	Conjugate:	Unconjugated
Calculated MW:	35kDa	Isotype:	IgG

Immunogen Information

Background: NIMA-related kinases share high amino acid sequence identity with the gene product of the *Aspergillus nidulans* 'never in mitosis A' gene, which controls initiation of mitosis.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: NIMA-related kinase 7; NEK7

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 203-302 of human NEK7 (Q8TDX7).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.