

CAB19820

Product Information

Product SKU:	CAB19820	Gene ID:	221656	Size:	20uL, 100uL
Clone No:	ARC2347	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	92kDa	Conjugate:	Unconjugated
Calculated MW:	92kDa	Isotype:	IgG

Immunogen Information

Background: Flavin-dependent histone demethylases, such as KDM1B, regulate histone lysine methylation, an epigenetic mark that regulates gene expression and chromatin function (Karytinis et al., 2009 [PubMed 19407342]).

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: AOF1; LSD2; C6orf193; LSD2/KDM1B/AOF1

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human LSD2/KDM1B/AOF1 (NP_694587.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.