

CAB19825

Product Information

Product SKU:	CAB19825	Gene ID:	259217	Size:	20uL, 100uL
Clone No:	ARC2351	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	75kDa	Conjugate:	Unconjugated
Calculated MW:	75kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable ATP binding activity. Located in extracellular exosome.

Recommended Dilution: WB,1:500 - 1:2000

Synonyms: HSPA12A ; HSPA12A

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 50-150 of human HSPA12A (O43301).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.