

CAB19827

Product Information

Product SKU:	CAB19827	Gene ID:	286887	Size:	20uL, 100uL
Clone No:	ARC2352	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	60kDa	Conjugate:	Unconjugated
Calculated MW:	60kDa	Isotype:	IgG

Immunogen Information

Background: Keratins are intermediate filament proteins responsible for the structural integrity of epithelial cells and are subdivided into epithelial keratins and hair keratins. The type II keratins are clustered in a region of chromosome 12q13.

Recommended Dilution: WB,1:500 - 1:2000 IHC-P,1:50 - 1:200

Synonyms: K6E; KRT6E; PPKNEFD; Cytokeratin 6 (KRT6)

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 465-564 of human Cytokeratin 6 (KRT6) (P48668).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,0.05% BSA,50% glycerol,pH7.3.