

CAB19971

---

## Product Information

<b>Product SKU:</b>	CAB19971	<b>Gene ID:</b>	8209	<b>Size:</b>	20uL, 100uL
<b>Clone No:</b>	-	<b>Host Species:</b>	Rabbit	<b>Reactivity:</b>	Human,Mouse,Rat

---

## Additional Information

<b>Observed MW:</b>	26kDa	<b>Conjugate:</b>	Unconjugated
<b>Calculated MW:</b>	24kDa/28kDa	<b>Isotype:</b>	IgG

---

## Immunogen Information

**Background:** This gene encodes a potential mitochondrial protein that is a member of the DJ-1/Pfpl gene family. This protein is overexpressed in fetal Down syndrome brain. Alternate splicing results in multiple transcript variants.

**Recommended Dilution:** WB,1:500 - 1:2000

**Synonyms:** ES1; HES1; KNPH; KNPI; GT335; GATD3A; C21orf33; 33

**Purification Method:** Affinity purification

**Immunogen:** Recombinant fusion protein containing a sequence corresponding to amino acids 42-268 of human C21orf33 (NP\_004640.4).

**Storage:** Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.