

CAB20147

Product Information

Product SKU: CAB20147	Gene ID: 10567	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 21kDa	Conjugate: Unconjugated
Calculated MW: 21kDa	Isotype: IgG

Immunogen Information

Background: Enables identical protein binding activity. Located in membrane.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: PRA1; YIP3; PRAF1; RABAC1

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-185 of human RABAC1 (NP_006414.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.