

CAB20148

Product Information

Product SKU:	CAB20148	Gene ID:	114785	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	101kDa	Conjugate:	Unconjugated
Calculated MW:	101kDa	Isotype:	IgG

Immunogen Information

Background: Enables chromatin binding activity. Located in chromocenter; fibrillar center; and nucleoplasm. Implicated in autism spectrum disorder.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 880-980 of human MBD6 (NP_443129.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.