

CAB20150

Product Information

Product SKU: CAB20150	Gene ID: 196513	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Rat

Additional Information

Observed MW: 75kDa	Conjugate: Unconjugated
Calculated MW: 68kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a member of a family of proteins that function in removing the 5' cap from mRNAs, which is a step in regulated mRNA decay. This protein localizes to cytoplasmic foci which are the site of mRNA breakdown and turnover. Alternative splicing results in multiple transcript variants.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: DCP1; DCP1B

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 386-617 of human DCP1B (NP_689853.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.