

CAB20151

Product Information

Product SKU: CAB20151	Gene ID: 199964	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: 24kDa	Conjugate: Unconjugated
Calculated MW: 22kDa	Isotype: IgG

Immunogen Information

Background: Predicted to be integral component of membrane.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 91-210 of human TMEM61 (NP_872338.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.