

CAB20155

Product Information

Product SKU:	CAB20155	Gene ID:	23067	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	250kDa	Conjugate:	Unconjugated
Calculated MW:	213kDa	Isotype:	IgG

Immunogen Information

Background: SET1B is a component of a histone methyltransferase complex that produces trimethylated histone H3 at Lys4 (Lee et al., 2007 [PubMed 17355966]).

Recommended Dilution: WB, 1:100 - 1:500

Synonyms: KMT2G; Set1B; IDDSELD; SETD1B

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1673-1770 of human SETD1B (NP_001340274.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.