

CAB20157

Product Information

Product SKU: CAB20157	Gene ID: 2846	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 42kDa	Conjugate: Unconjugated
Calculated MW: 42kDa	Isotype: IgG

Immunogen Information

Background: This gene encodes a member of the lysophosphatidic acid receptor family. It may also be related to the P2Y receptors, a family of receptors that bind purine and pyrimidine nucleotides and are coupled to G proteins. The encoded protein may play a role in monocytic differentiation.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: LPA4; P2Y9; GPR23; P2RY9; P2Y5-LIKE; LPAR4

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 311-370 of human LPAR4 (Q99677).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.