

CAB20166

Product Information

Product SKU: CAB20166	Gene ID: 55106	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human

Additional Information

Observed MW: 67kDa	Conjugate: Unconjugated
Calculated MW: 67kDa	Isotype: IgG

Immunogen Information

Background: Predicted to act upstream of or within negative regulation of cell population proliferation.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: SLFN3; SLFN12

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 320-390 of human SLFN12 (NP_001275938.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.