

CAB20168

Product Information

Product SKU:	CAB20168	Gene ID:	55192	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	35kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to enable RNA binding activity. Predicted to act upstream of or within negative regulation of transcription by RNA polymerase II and toxin transport. Predicted to be located in cytoplasm and nucleus.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: -

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-304 of human DNAJC17 (NP_060633.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.