

CAB20173

Product Information

Product SKU:	CAB20173	Gene ID:	64855	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	95kDa	Conjugate:	Unconjugated
Calculated MW:	84kDa	Isotype:	IgG

Immunogen Information

Background: Enables transcription coactivator activity. Involved in several processes, including gonadotropin secretion; positive regulation of transcription regulatory region DNA binding activity; and regulation of cellular macromolecule biosynthetic process. Located in several cellular components, including adherens junction; cytosol; and nucleoplasm.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: OC58; MEG-3; C9orf88; FAM129B; MINERVA; bA356B19.6

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 447-746 of human FAM129B (NP_073744.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.