

CAB20180

Product Information

Product SKU:	CAB20180	Gene ID:	83903	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	100kDa	Conjugate:	Unconjugated
Calculated MW:	88kDa	Isotype:	IgG

Immunogen Information

Background: Enables ATP binding activity and histone kinase activity (H3-T3 specific). Involved in histone H3-T3 phosphorylation involved in chromosome passenger complex localization to kinetochore; intracellular signal transduction; and mitotic sister chromatid cohesion. Located in centrosome; nucleoplasm; and spindle.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: GSG2; HASPIN

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 290-470 of human HASPIN (NP_114171.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.