

CAB20182

Product Information

Product SKU: CAB20182	Gene ID: 84549	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 34kDa	Conjugate: Unconjugated
Calculated MW: 35kDa	Isotype: IgG

Immunogen Information

Background: Enables RNA binding activity. Predicted to be involved in maturation of 5.8S rRNA and maturation of LSU-rRNA. Located in nucleolus.

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: RBM13; MAK16L; MAK16

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 30-120 of human MAK16 (NP_115898.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.