

CAB20184

Product Information

Product SKU:	CAB20184	Gene ID:	8290/8350	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat,Other (Wide Range Predicted)

Additional Information

Observed MW:	17kDa	Conjugate:	Unconjugated
Calculated MW:	15kDa	Isotype:	IgG

Immunogen Information

Background: Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. This structure consists of approximately 146 bp of DNA wrapped around a nucleosome, an octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H3 family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. This gene is found in the large histone gene cluster on chromosome 6p22-p21.3.

Recommended Dilution: DB,1:500 - 1:1000 WB,1:500 - 1:1000 ChIP,5µg antibody for 5µg-10µg of Chromatin

Synonyms: H3/A; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FA; H3C10; H3C11; H3C12; HIST1H3A; Acetyl-Histone H3-K27

Purification Method: Affinity purification

Immunogen: A synthetic acetylated peptide around K27 of human Histone H3 (NP_003520.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.