

CAB20191

Product Information

Product SKU:	CAB20191	Gene ID:	11433347	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Medicago truncatula

Additional Information

Observed MW:	15kDa	Conjugate:	Unconjugated
Calculated MW:	15kDa	Isotype:	IgG

Immunogen Information

Background: Provides oxygen to the bacteroids. This role is essential for symbiotic nitrogen fixation.

Recommended Dilution: WB, 1:1000 - 1:2000

Synonyms: LB2; MTLb2; MtLBs

Purification Method: Affinity purification

Immunogen: Two synthetic peptides corresponding to the sequence within amino acids 7-21 and 100-113 of Medicago truncatula Leghemoglobin 2 (XP_003588689.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.