

CAB20212

Product Information

Product SKU: CAB20212	Gene ID: 220988	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 40kDa	Conjugate: Unconjugated
Calculated MW: 40kDa	Isotype: IgG

Immunogen Information

Background: Enables RNA binding activity. Predicted to be involved in mRNA splicing, via spliceosome. Located in nucleoplasm. Part of catalytic step 2 spliceosome.

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: FBRNP; HNRPA3; D10S102; 2610510D13Rik; HNRNPA3

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 100-200 of human HNRNPA3 (NP_001317178.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Contact Details | Dublin, Ireland

Email: techsupport@assaygenie.com | **Web:** www.assaygenie.com

Copyright © 2020 Reagent Genie, All Rights Reserved. All information / detail is correct at time of going to print.