

CAB20216

Product Information

Product SKU:	CAB20216	Gene ID:	3383	Size:	100uL
Clone No:	ARC5129-01	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	Refer to figures	Conjugate:	Unconjugated
Calculated MW:	58kDa	Isotype:	IgG

Immunogen Information

Background: This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor.

Recommended Dilution: Elisa, 1:500 - 1:2000 FC, 1:50 - 1:200 IF/ICC, 1:50 - 1:200

Synonyms: BB2; CD54; P3.58; ICAM-1/CD54

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 28-480 of human ICAM1.(NP_000192.2)

Storage: Store at 4°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% thimerosal, pH7.3.