

CAB20222

Product Information

Product SKU:	CAB20222	Gene ID:	1045	Size:	20uL, 100uL
Clone No:	ARC3070	Host Species:	Rabbit	Reactivity:	Human,Mouse

Additional Information

Observed MW:	38kDa	Conjugate:	Unconjugated
Calculated MW:	34kDa	Isotype:	IgG

Immunogen Information

Background: This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded protein is a major regulator of intestine-specific genes involved in cell growth and differentiation. This protein also plays a role in early embryonic development of the intestinal tract. Aberrant expression of this gene is associated with intestinal inflammation and tumorigenesis.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: CDX3; CDX-3; CDX2/AS; X2

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CDX2 (Q99626).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.