

CAB20245

Product Information

Product SKU: CAB20245	Gene ID: 108161	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Human,Mouse,Rat

Additional Information

Observed MW: 39kDa	Conjugate: Unconjugated
Calculated MW: 40kDa	Isotype: IgG

Immunogen Information

Background: Predicted to be involved in chromatin organization. Predicted to be located in intercellular bridge; midbody; and nucleoplasm. Predicted to be active in nucleus. Is expressed in primary spermatocyte; seminiferous tubule; and spermatid. Orthologous to human FAM50B (family with sequence similarity 50 member B).

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: X5L; D0H6S2654E; XAP-5-like; Fam50b

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 70-150 of mouse Fam50b (NP_620085.2).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.