

CAB20275

Product Information

Product SKU:	CAB20275	Gene ID:	241556	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Mouse,Rat

Additional Information

Observed MW:	28kDa	Conjugate:	Unconjugated
Calculated MW:	28kDa	Isotype:	IgG

Immunogen Information

Background: Predicted to be located in membrane. Predicted to be integral component of membrane. Predicted to be integral component of plasma membrane. Is expressed in several structures, including alimentary system; brain; genitourinary system; immune system; and integumental system. Orthologous to human TSPAN18 (tetraspanin 18).

Recommended Dilution: WB,1:500 - 1:1000

Synonyms: 6720430O15; 2610042G18Rik; Tspan18

Purification Method: Affinity purification

Immunogen: A synthetic peptide corresponding to a sequence within amino acids 149-248 of mouse Tspan18 (NP_899003.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.