

CAB20327

Product Information

Product SKU:	CAB20327	Gene ID:	107276161	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Oryza sativa

Additional Information

Observed MW:	27kDa	Conjugate:	Unconjugated
Calculated MW:	27kDa	Isotype:	IgG

Immunogen Information

Background: Probable transcription factor involved in the regulation of flowering time under long day (LD conditions). Plays a major role as repressor of flowering. Controls flowering time by negatively regulating the expression of EHD1 and HD3A.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: Ghd7; GHD7

Purification Method: Affinity purification

Immunogen: Recombinant protein of Oryza sativa GHD7.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.