

CAB20329

Product Information

Product SKU: CAB20329	Gene ID: 107278653	Size: 20uL, 100uL
Clone No: -	Host Species: Rabbit	Reactivity: Oryza sativa

Additional Information

Observed MW: 50kDa	Conjugate: Unconjugated
Calculated MW: 46kDa	Isotype: IgG

Immunogen Information

Background: Putative transcription regulator that controls rice tillering by initiating axillary buds and promoting their outgrowth. Rice tiller is a specialized grain-bearing branch that is formed on the unelongated basal internode and grows independently of the mother stem (culm) by means of its own adventitious roots.

Recommended Dilution: WB, 1:500 - 1:2000

Synonyms: MOC1

Purification Method: Affinity purification

Immunogen: A synthetic peptide of Oryza sativa MOC1.

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.