

CAB20335

Product Information

Product SKU:	CAB20335	Gene ID:	2267	Size:	20uL, 100uL
Clone No:	ARC5134-01	Host Species:	Rabbit	Reactivity:	Human

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	36kDa	Isotype:	IgG

Immunogen Information

Background: Fibrinogen-like 1 is a member of the fibrinogen family. This protein is homologous to the carboxy terminus of the fibrinogen beta- and gamma- subunits which contains the four conserved cysteines of fibrinogens and fibrinogen related proteins. However, this protein lacks the platelet-binding site, cross-linking region and a thrombin-sensitive site which are necessary for fibrin clot formation. This protein may play a role in the development of hepatocellular carcinomas. Four alternatively spliced transcript variants encoding the same protein exist for this gene.

Recommended Dilution: WB, 1:500 - 1:1000

Synonyms: HPS; HFREP1; HP-041; LFIRE1; LFIRE-1; FGL1

Purification Method: Affinity purification

Immunogen: Recombinant protein of human FGL-1 Full length (NP_004458.3).

Storage: Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 0.05% BSA, 50% glycerol, pH 7.3.