

CAB20338

Product Information

Product SKU:	CAB20338	Gene ID:	69806	Size:	20uL, 100uL
Clone No:	-	Host Species:	Rabbit	Reactivity:	Human,Mouse,Rat

Additional Information

Observed MW:	35kDa	Conjugate:	Unconjugated
Calculated MW:	35kDa	Isotype:	IgG

Immunogen Information

Background: Enables zinc ion transmembrane transporter activity. Predicted to be involved in zinc ion transmembrane transport. Predicted to act upstream of or within zinc ion transport. Located in Golgi apparatus; nucleus; and plasma membrane. Orthologous to human SLC39A11 (solute carrier family 39 member 11).

Recommended Dilution: WB,1:500 - 1:1000 IHC-P,1:50 - 1:200

Synonyms: ZIP-11; 1810074D23Rik; Slc39a11

Purification Method: Affinity purification

Immunogen: Recombinant fusion protein containing a sequence corresponding to amino acids 101-174 of mouse Slc39a11 (NP_001159975.1).

Storage: Store at -20°C. Avoid freeze / thaw cycles.Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.